

Policy & Procedure (P&P)

Policy Title :

Labeling of Blood Components Units

| Department | Index No. | Scope |
|-------------------------|----------------------|----------------------|
| Laboratory & Blood Bank | LAB-068 | All Blood Bank Staff |
| Issue Date | Revision NO | Effective Date |
| 1439/9/8 | 2 | 1440/07/23 |
| Review Due Date | Related Standard NO. | Page Number# |
| 1442/07/23 | CBAHI (LB. 53) | Page 1 of 5 |

01. Policy:

- 01.1. The blood bank develops an initial labeling on each blood component unit after completion of the donor testing and not before the discard of unacceptable units.

02. Definition :

- 01.1. N/A

03. Purpose :

- 01.1. The objective of the blood component labeling is to reduce the danger of incompatible transfusions caused through human errors by presenting important information in a clear, logical and easily recognizable format.

04. Procedure :

04.1. Material Required:

- 04.1.1. The label model on Microsoft publisher
- 04.1.2. ZEBRA printer
- 04.1.3. Adhesive labels for all components

04.2. Procedure

- 04.2.1. After collection and processing whole blood and blood components, all the units remain in quarantine storage areas (Unscreened blood bank refrigerator, deep freezer and platelet incubator and agitator).
- 04.2.2. Once all the reports of blood group and TTI testing are ready, the blood bank staff place the bags

on a table in chronological order.

- 04.2.3. Segregate those which are found reactive for any TTI or found unsuitable for use and keep them in the area for disposal. Leave those found suitable for use on the bench for labeling.
- 04.2.4. On the computer choose the specific label model on Microsoft publisher for the specific blood component (PRBCs or Platelets or FFP or Cryoprecipitate) Write carefully the unit number, date of collection and expiry date and the Blood group and the phenotype on each label as per the grouping register records.
- 04.2.5. Date of collection and date of expiry is very important. The expiry date depends on the type of bag and component.

The day of blood collection is considered the day zero for calculating the expiry dates.


04.3. The initial labeling:

- 04.3.1. The initial labeling requirement include:
- 04.3.2. The identification of the collection facility: we are using the Letter Q to refer to Qunfudah General Hospital Blood Bank followed by the year (for example Q39) and the logo of the Qunfudah General Hospital.
- 04.3.3. The product names.
- 04.3.4. The unit number: for example, Q39-1520
- 04.3.5. The ABO-Rh
- 04.3.6. The collection date and time
- 04.3.7. The expiration date and time
- 04.3.8. The mention that the unit is non-reactive for TTD (HBSAg, anti HBC, antiHCV, HIVP24, anti HIV1,2, anti HTLV1,2, VDRL, NAT and malaria)
- 04.3.9. The mention that the antibody screening is negative.

04.4. Initial labeling for PRBCs:

The blood bank staff will add:


- The product name: PRBCs
- The phenotype Rh/Kell
- The mention that it is Leukodepleted
- If the unit is irradiated then the mention IRRADIATED will be added
- The blood bank staff will add the mention SICKLING TEST NEGATIVE for the unit prepared to the SCD patients.
- This is a sample for our labeling sticker on the PRBCs units:

| | | |
|---|---|--|
| PRBCs- Unit No: Q39-1548 | |  |
| Collection Date 20-5-2018 | B+ Rh Positive CcDee Kneg LEUKODEPLETED Additive Solution: S.AGM:100 ml Stored at -1 C to -6 C | |
| Expiry Date 1-7-2018 @12 am | | |
| NON REACTIVE for : HBsAg, anti HBC, anti HCV, HIV P24, anti HIV 1,2, anti HTLV 1,2, VDRL , NAT and malaria Antibody screening Neg | | |

04.4.1. **Initial labeling for the platelet's units:**

The initial labeling requirement include:


- The identification of the collection facility: we are using the Letter Q to refer to Qunfudah General Hospital Blood Bank followed by the year (for example Q39) and the logo of the Qunfudah General Hospital.
- The product name: Platelets concentrate or Platelets Apheresis
- The mention that it is Leukodepleted if it is leukodepleted
- The mention that Bacterial contamination was not detected at the date and time of issue
- This is a sample for our labeling sticker on the Platelets units:

| | | |
|--|--|--|
| Platelets- Unit No: Q39-1549 | |  |
| Collection Date 20-5-2018 | O+ Rh positive LEUKODEPLETED stored at -20-24 C Bacterial contamination not detected at the date and time of issue | |
| Expiry Date 25-5-2018 @12am | | |
| NON REACTIVE for : HBsAg, anti HBC, anti HCV, HIV P24, anti HIV 1,2, anti HTLV 1,2, HBV DNA, HCV RNA, HIV RNA VDRL and malaria | | |

04.4.2. **Initial labeling for the fresh frozen plasma units:**


- The product name: FFP (fresh frozen plasma)

- This is a sample for our labeling sticker on the FFP units:

| | | |
|---|--------------------|--|
| FFP- Unit No: Q39-1549 | |  |
| Collection Date 20-5-2018 | O+ | |
| Expiry Date 20-5-2019 @12am | Rh positive | |
| NON REACTIVE for : HBsAg, anti HBC, antiHCV, HTP24, anti HIV1.2, anti HTLV1.2, HBV DNA, HCV RNA, HIV RNA VDRL and malaria | | stored at—15 C |

04.4.3. **Initial labeling for the cryoprecipitate units:**

- The product name: cryoprecipitate
- This is a sample for our labeling sticker on the cryoprecipitate units:

| | | |
|---|--------------------|--|
| Cryo - Unit No: Q39-1549 | |  |
| Collection Date 20-5-2018 | A+ | |
| Expiry Date 20-5-2019 @12am | Rh positive | |
| NON REACTIVE for : HBsAg, anti HBC, antiHCV, HTP24, anti HIV1.2, anti HTLV1.2, HBV DNA, HCV RNA, HIV RNA VDRL and malaria | | stored at—15 C |

- 04.5. After the bags are labeled, a second technician double checks the number and blood groups on the bags tallying them with the records.

05. Responsibilities :

- 05.1. All Blood Bank Staff of Al-Qunfudah General Hospital.



وزارة الصحة

Ministry of Health

مستشفى القنفذة العلم

06. Equipment & Forms

06.1. N/A

07. Attachment :

07.1. N/A

08. Reference

08.1. The Technical manual of the American Association of Blood Banks.

Preparation , Reviewing & Approval Box

| | NAME | POSITION | SIGN & STAMP | DATE |
|----------------------|-----------------------|--------------------|--------------|-----------|
| Prepared By | Dr RAJA NACER SASSI | Head of Blood Bank | | |
| Reviewed By | Mr. ABDULHADI ASHIRI | Lab & B.Bank HOD | | |
| Document Reviewed By | Ms. SADIAH ALMAHMOUDI | TQM Director | | 122/9/5 |
| Reviewed By | Dr. AGEEL ALGANIMI | Medical Director | | |
| Approved By | Dr. ABDULLAH ALJABRI | Hospital Director | | 1430/7/14 |



COLLECTION DAY -DAY 0-

COLLECTION DATE.....

| PERFORMED AND ACCEPTABLE? | YES | NO | NIU |
|--|-----|----|-----|
| RECEPTION AND DONATION SECTION | | | |
| HEMOCUE QUALITY CONTROL | | | |
| BLOOD SCHAKER QUALITY CONTROL | | | |
| DONOR HISTORY QUESTIONNAIRES FILLED AND REVIEWED | | | |
| CRASH CART DAILY CHECK | | | |
| TRIMA DAILY MAINTENANCE | | | |

COMMENTS:.....

.....

.....

I certify that I have reviewed all above records, listed all deficiencies and notified the supervisor.

Signature reception technician

signature donation technician

COLLECTION DAY -DAY 0-

COLLECTION DATE.....Components Separation date.....

| PERFORMED AND ACCEPTABLE? | YES | NO | NIU |
|-------------------------------------|-----|----|-----|
| REVEOS DAILY CHECK | | | |
| REFRIGERATED CENTRIFUGE DAILY CHECK | | | |
| WELDER DAILY CHECK | | | |
| SEALER DAILY CHECK | | | |
| BALANCE DAILY CHECK | | | |
| ORTHOVISION DAILY QC | | | |
| TANGO INFINITY DAILY QC | | | |
| ID GEL CARD DAILY QC | | | |

COMMENTS:.....

.....

.....

I certify that I have reviewed all above records, listed all deficiencies and notified the supervisor.

Signature separation technician

signature X-match technician

INVESTIGATIONS DAYS – DAYS 1/2-

Date:.....

| PERFORMED AND ACCEPTABLE? | YES | NO | NIU |
|---|-----|----|-----|
| DONOR ABO-RH GROUPING + PHENOTYPE+ ANTIBODY SCREENING RESULTS REVIEWED BY SENIOR | | | |
| DONOR BLOOD GROUPING FROM SEGMENTS REVIEWED BY SENIOR AND MATCHED WITH THOSE FROM SAMPLES | | | |
| DONOR MALARIA RESULTS AUTHORIZED BY SENIOR | | | |
| DONOR SEROLOGY + NAT RESULTS APPROVED AND RECEIVED FROM SEROLOGY DEPARTMENT | | | |
| OTHER TESTS : Du AND SICKLING TESTS | | | |
| PLATELETS BACTERIAL DETECTION | | | |
| DAILY VISUAL INSPECTION | | | |

COMMENTS:.....

.....

I certify that I have reviewed all above records and are all acceptable .

I authorize the labeling and release of the blood components from this collection date.

Signature blood bank consultant

Signature blood bank specialist

LABELS STICKERS

SERIAL NUMBERS FROM

TO

RECORDS REVIEWED BY

LABELS STICKERS PRINTED BY

DOUBLE CHECKED BY

DATE